



How to use your

Rēd Comb

LED light care for the scalp

Before you start

- Hair should be dry and detangled.
- Divide the hair into approximately 4–6 sections for easier and more even use.

Combing technique

1. SECTION

Create a clean part and move the remaining hair aside so the scalp is easily accessible.

2. POSITION

Place the comb teeth directly against the scalp.

Guide the device so that the LED light sources remain close to the scalp during use.

3. GLIDE

Move the comb slowly and smoothly across the selected area. Work systematically to provide even light coverage across the section.

4. REPEAT

Move to the next section and continue until the desired scalp area has been covered.

Session duration

- Approximately 2–3 minutes per section
- Total session approximately 10–15 minutes
- Follow the recommended duration and frequency stated in the user manual

Select mode

RED: Visible red LED light

SECOND LIGHT MODE: According to the specification of the individual device

COMBO: Simultaneous use of the available light modes, where supported by the device

Frequency

- Incorporate Rēd Comb into your regular hair-care and wellness routine.
- Use according to the recommended frequency stated in the user manual
- Follow the recommended duration for each session
- Do not exceed the manufacturer's instructions for frequency or duration of use

What to expect

Rēd Comb provides a simple way to incorporate LED light care into a regular routine focused on the scalp.

Regular use can help you maintain a consistent hair-care routine and give different areas of the scalp equal attention.

Suitable for

- Regular LED light care for the scalp
- Home hair-care and wellness routines
- Convenient use on selected areas of the scalp
- Salons, spas and professional wellness environments

Important note

- Rēd Comb is intended for cosmetic, wellness and personal light-care use.
- Use the device only as described in the user manual.
- Do not use the device on irritated, injured or otherwise problematic skin without appropriate professional advice.
- Stop using the device if you experience an unexpected or uncomfortable reaction.

Tips

- Use on dry, detangled hair.
- Dividing the hair into sections helps provide more even use.
- Move the device slowly and without excessive pressure.
- After use, continue with your usual hair-care routine.
- Clean the device after use according to the instructions in the user manual.